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Milk bank pasteurizer

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Wholly made of 18/10 AISI 304 stainless steel. Upper smooth surface and upper level at the back. One heating and cooling sink, sizes as per drawing, insulated. Upper hinged double-wall cover complete with a gas spring keeping it under safe conditions when it is open. Hydraulic plant made up of solenoid valves and water charge duct \varnothing 3/4", solenoid valves and anti-steam pipe coil on the edge of the sink \varnothing 1/2". Lower motorised drain with motorised ball valve \varnothing 1". Electric installation made up of a main switch, control panel complete with a backlit liquid crystal display to display the phases of the operation cycle, control keys and mushroom-shaped button for reset or emergency. The electronic board is arranged inside the lower compartment to receive/distribute input/output signals, a safety thermoregulator, remote switches to control heating and cooling resistances (3 2500 W resistances). End-of-cycle buzzer. General disconnecter with interlock of the electric panel. Refrigerating compressor unit installed on the lower wall, of a water-cooled type, 2 Hp power. Chiller refrigerating plant with water recirculation cooled in a temperate sink of use + 2 +4°C. Type gas R134A.

SAMPLE TEMPERATURE SENSOR:

The sample temperature sensor is intended to sample the pasteurizing cycle. It is designed for insertion in a control bottle. This sensor is operating in place of the sensor in the sink. It is activated by inserting it into the plug on the control display and it is intended to control the work cycle. The temperature values measured by the sensor during the operation are recorded and stored by the controller. Afterwards it will be possible to check every single work cycle by means of a management software.

Work cycle:

- hot water charge in automatic mode, regulated by a level sensor
- heating up to the pre-set temperature (62-64°C) for the pasteurizing time established (30 minutes)
- hot water discharge
- sink filling in automatic mode, regulated by a level sensor
- water and nursing bottle cooling up to the pre-set temperature (minimum admitted value +2°C.)
- sound / light end-of-cycle warning signal

Lower piece of furniture closed on either side, opened at the back. Assembly on adjustable feet.

European directives applied: 89/392 EEC (machines), 73/23 EEC (low voltage), 89/336 EEC (electromagnetic compatibility). European standards of reference EN 60601-1.

Scotch brite finish.

SOFTWARE INTENDED TO MANAGE AND DOWNLOAD THE DATA FOR THE PASTEURIZER:

The "Pasteur ver. 1.0.4" software realised to record, download, display and print the cycle data for the Metalarredinox pasteurizer is operating on a Windows™ platform. It is activated by a hardware key connected to a parallel or a USB port of the PC.

The pasteurizer during the work cycle is recording the temperature measured by the active temperature sensor at regular intervals (every 5 min.) (sensor of a fixed type immersed into the water bath or mobile sensor inserted into the sample bottle).

By connecting a personal computer with the pasteurizer at the end of a work cycle (the connection occurs between a serial plug of the PC and the special plug on the display), you can download the series of data acquired.

From the PC you can select the specific pasteurizing cycle (you can recognise it by the start date and time and by a progressive number), display and print: the temperature diagram during the cycle, the series of data in the form of a table (time, temperature and cycle phase, heating of the pasteurizing process or cooling) and the setup of machine parameters.

PRINTER FOR HACCP PASTEURIZER





Chemical paper printer intended to record the pasteurizing and cooling cycle, not depending upon the pasteurizer. It is equipped with a temperature sensor of the PT 100 type to be inserted into the sample bottle by means of a passage complete with a gasket. The printer is made up of a stainless steel case to arrange it easily on a working surface. It uses 57.5 mm. rolls of chemical paper. The sampling interval is set to 5 minutes (you can change the reading interval in the setup menu).

The output ticket will indicate not only the start-of-record and end-of-record date and time and the measured data but also supply the following information: equipment, product, operator's name, approved.

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METALARREDINOX consistently with its policy of continuous improvement on its own products, reserves the right to make changes at any time without prior notice. Therefore, all related data must be regarded as approximate and unbinding.